




PATIENT RESULTS
0123456-78
Patient, John A.

701 Congressional Blvd., Suite 360
Carmel, IN 46032
Phone: (800) 571-1948
Fax: (317) 573-263
CLIA ID: 0000-ABC-1234

Client: Orchard Office	Patient Name: Patient, John A.	Requisition #: 1234-56
Provider: Orchard Orchard	Date of Birth: 2/28/1990	Collected: 3/24/2015
Phone: (800) 571-5835	Patient ID: 000-2345-12	Received: 3/24/2015
Fax: (317) 573-2633	Sample ID: 0123456-78	Approval Date: 3/24/2015

Confirmed Analytes CLON, MORPH, TEMAZ	COD	Prescribed Medication(s) Kadian (Morphine), Lortab/Norco (Hydrocodone)
---	-----	--

Summary Of All Laboratory Testing						
Analyte(s)	Physician POCT	Lab EIA	LC/MS/MS Result	Level (ng/mL)	Cutoff (ng/mL)	Comments
Opiates						Significant historical finding of multiple mutations on CYP450 3A4/3A5 and designated as having poor metabolism of opiate medications.
Codeine		Pos	Confirmed	82	50 ng/mL	
6-MAM	Inconsistent	Neg	Confirmed	97	50 ng/mL	
Morphine		Pos	Negative	<50	50 ng/mL	
Hydrocodone		Neg	Negative	<50	50 ng/mL	
Hydromorphone	Inconsistent	Neg	Negative	<50	50 ng/mL	



Benzodiazepines						
Analyte(s)	Physician POCT	Lab EIA	LC/MS/MS Result	Level (ng/mL)	Cutoff (ng/mL)	Comments
Clonazepam		N/A	Confirmed	25	20 ng/mL	
Alprazolam		N/A	Negative	<50	50 ng/mL	
Lorazepam		N/A	Negative	<20	20 ng/mL	
Nordiazepam		N/A	Negative	<20	20 ng/mL	
Oxazepam		N/A	Negative	<15	15 ng/mL	
Temazepam		N/A	Confirmed	69	20 ng/mL	



PATIENT RESULTS
0123456-78
Patient, John A.

701 Congressional Blvd., Suite 360
Carmel, IN 46032
Phone: (800) 571-1948
Fax: (317) 573-263
CLIA ID: 0000-ABC-1234

Client: Orchard Office	Patient Name: Patient, John A.	Requisition #: 1234-56
Provider: Orchard Orchard	Date of Birth: 2/28/1990	Collected: 3/24/2015
Phone: (800) 571-5835	Patient ID: 000-2345-12	Received: 3/24/2015
Fax: (317) 573-2633	Sample ID: 0123456-78	Approval Date: 3/24/2015

<u>Confirmed Analytes</u> CLON, MORPH, TEMAZ	COD	<u>Prescribed Medication(s)</u> Kadian (Morphine),Lortab/Norco (Hydrocodone)
---	-----	---

Summary Of All Laboratory Testing						
Analyte(s)	Physician POCT	Lab EIA	LC/MS/MS Result	Level (ng/mL)	Cutoff (ng/mL)	Comments
Opiates						Significant historical finding of multiple mutations on CYP450 3A4/3A5 and designated as having poor metabolism of opiate medications.
Codeine	Inconsistent	Pos	Confirmed	61	50 ng/mL	
6-MAM		Neg	Negative	<50	50 ng/mL	
Morphine		Pos	CONFIRMED	97	50 ng/mL	
Hydrocodone	Inconsistent	Neg	Negative	<50	50 ng/mL	
Hydromorphone		Neg	Negative	<50	50 ng/mL	
Benzodiazepines						
Clonazepam		N/A	Confirmed	25	20 ng/mL	
Alprazolam		N/A	Negative	<50	50 ng/mL	
Lorazepam		N/A	Negative	<20	20 ng/mL	
Nordiazepam		N/A	Negative	<20	20 ng/mL	
Oxazepam		N/A	Negative	<15	15 ng/mL	
Temazepam		N/A	Confirmed	69	20 ng/mL	



LABORATORY REPORT

701 Congressional Blvd.
Carmel, IN 46032
Phone: (800) 571-5835
Fax: (317) 573-2633

LOCATION DETAILS

Name: Orchard Office
 Account:
 Address: 701 Congressional Blvd.
 Suite 360
 Fax: (317) 573-2633
 Provider: Orchy Orchardson

PATIENT DETAILS

Name: Patient, John A.
 DOB: 5/23/1981
 Gender: M
 ID: 1

SPECIMEN DETAILS

Requisition Number:
 Lab Accession Number: 1047
 Date Collected: 1/30/2014
 Date Received by Lab: 1/30/2014 2:38 PM
 Date Reported: 00/00/00

Analyte Tested	Lab Result (Qualitative)	Lab Result (ng/mL)	Assay Cutoff (ng/mL)	Medication Comparison	Previous Result (T-1)	Detection Window
Benzodiazepines						
Clonazepam	Negative	0	<20			
Alprazolam	Negative	0	<50			
Lorazepam	Negative	0	<10			
Nordiazepam	Negative	0	>20			
Oxazepam	Negative	0	<15			
Temazepam	Negative	0	<20			
Opiates						
Codeine	Confirmed	82	50		Negative	
6-MAM	Confirmed	97	50	Inconsistent	5.00	This is a test comment.
Morphine	Negative	<50	50		Negative	
Hydrocodone	Negative	<50	50		Negative	
Hydromorphone	Negative	48	50	Inconsistent	Negative	
RESULTS EXPLANATION	This is an order choice comment. This section will continue throughout the length of the order choice section. It will extend the order choice section beyond the last test in the order choice if there is not enough room. Otherwise, the order choice section will end after the last test in the order choice.					



LABORATORY REPORT

701 Congressional Blvd.
Carmel, IN 46032
Phone: (800) 571-5835
Fax: (317) 573-2633

LOCATION DETAILS

Name: Orchard Office
Account:
Address: 701 Congressional Blvd.
Suite 360
Fax: (317) 573-2633
Provider: Orchy Orchardson

PATIENT DETAILS

Name: Patient, Jane A.
DOB: 4/4/1984
Gender: F
ID: D4444444

SPECIMEN DETAILS

Requisition Number:
Lab Accession Number: 1060
Date Collected: 5/12/2014
Date Received by Lab: 5/12/2014 1:45 PM
Date Reported: 00/00/00

CONSISTENT POSITIVE RESULTS - Suggests that the patient has been taking medication as prescribed

PRESCRIBED MEDICATION	EXPECTED METABOLITE	LEVEL	DETECTION WINDOW	COMMENTS
Duragesic	Fentanyl (MS)	115 ng/mL	24 hours (single intravenous dose). Higher doses & certain pathologies may extend the detection window.	The specimen tested positive for Fentanyl, Norfentanyl, which is consistent with patient having taken DURAGESIC.
	Norfentanyl (MS)	94 ng/mL	72 hours (single intravenous dose). Higher doses & certain pathologies may extend the detection window.	

INCONSISTENT NEGATIVE RESULTS - Suggests that the patient has not been taking prescribed medications

PRESCRIBED MEDICATION	EXPECTED METABOLITE	DETECTION WINDOW	COMMENTS
Xanax	Alprazolam	2 - 5 days after last dose. Higher doses & certain pathologies may extend the detection window.	The specimen tested negative for Alprazolam, which suggests the patient is not currently taking XANAX (within the detection window)

INCONSISTENT POSITIVE RESULTS - Suggests the patient may be taking illicit or unprescribed substances

DETECTED SUBSTANCE	LEVEL	DETECTION WINDOW	COMMENTS
Noroxycodone	147 ng/mL	Typically follow oxycodone concentrations	
Oxycodone	184 ng/mL	1 - 2 days (immediate release); 2 - 4 days (extended release)	
Oxymorphone	119 ng/mL	2 - 3 days (immediate release); 2 - 4 days (extended release)	

SPECIMEN VALIDITY - Tests to determine if the sample has been altered or results obscured

TEST	RESULT	LEVEL	REFERENCE RANGE	UNITS	COMMENTS
Urine Creatinine	Normal	25	>20	mg/dL	
Oxidants	Normal	178	>200	ug/mL	
Urine pH	Normal	8.1	4.5-9.5		
Urine Sp. Gravity	Normal	1.039	1.003-1.040		

This text is part of the footer. In addition to the end of report text, you can attach text to the footer displayed at the bottom of every page. The footer can be as long or as short as you need it to be, and can contain various demographic information about the order as well as block text.



700 Congressional Blvd.
Suite 360
Carmel, IN 46032

Orchard Office
701 Congressional Blvd.
Suite 360
Orchy Orchardson

Collected Date: 11/6/2013
Analysis Requested: InconRes, ConRes,

Tox Number: 0003
Reported Date: 11/7/2013
Specimen ID: 1043

INCONSISTENT RESULTS for urine collected 11/6/2013			
DRUG	LEVEL*	INTERPRETATION AND MEDICATION INFORMATION	
Hydromorphone	<20	Quantity is less than expected based on dosage regimen.	

CONSISTENT RESULTS for urine collected 11/6/2013			
DRUG	LEVEL*	INTERPRETATION AND MEDICATION INFORMATION	
Hydrocodone	397	Result is expected based on prescription information.	
Dihydrocodeine	210	Result is expected based on parent drug information.	

VALIDITY TESTING for urine collected 11/6/2013			
SPECIMEN VALIDITY	RESULT	NORMAL RANGE	INTERPRETATION
Urine Creatinine	41	>20	
Urine pH	6.9	3.0-11.0	
Urine Sp. Gravity	1.020	1.003-1.035	

Certifying Scientist: Orchy C. Orchardson **For more information, call (800) 571-5835.**
 * All levels have been normalized to creatinine (ng/mg)

HISTORICAL RESULTS for urine collected			
DRUG	LEVEL*	RESULT	INTERPRETATION
Hydrocodone	Negative		
Hydromorphone	Negative		
Dihydrocodeine	Negative		
Pregabalin	Negative		

OrchardSoftware



Harvest the Power

Sample reports generated from Orchard Software test database.

Note that these are sample reports generated with tools available in Orchard Software's laboratory solutions. Specific report formatting assistance from Orchard Software may require a professional services agreement.

All patient and provider names and data contained within these sample reports are fictitious and do not represent real names or data.

© Copyright 2015 Orchard Software Corporation. All rights reserved.